Lab 10: Earth's Magnetic Field

Philip Kim

April 14, 2021

- 1. Record the initial dip angle $\theta_0=36^\circ$
- 2. Set source to 4V.

| Table 1: High-Pass Filter | | | | | | |
|---------------------------|------------|------------|------------|-------------|-------------|-------------|
| Resistance | 20Ω | 40Ω | 75Ω | 150Ω | 180Ω | 200Ω |
| Current i | 0.122A | 0.0733A | 0.0442A | 21.24 mA | 0.02124A | 0.01824A |
| Dip Angle Θ_i | -71° | -49° | -14° | 6° | 12° | 15° |
| Calculated B_i | 1 | 2 | 3 | 4 | 5 | 6 |

- 3. Record the Helmholtz coil radius: R = 9.75cm \rightarrow 0.0975m
- 4. Record the Helmholtz coil number of turns: N=128
- 5. Calculations: $(B_i = \frac{8N_{\mu_0}I_i}{R\sqrt{125}}, where \mu_0 = 4\pi \times 10^{-7}Tm/A)$
 - Plot $tan\theta_i$ vs B_i with straight line. Deduce the values of B_V and B_H from the graph.

$$B_V =$$

$$B_H =$$

Calculate B_E

Lookup value of \mathcal{B}_{E}